



PRODUCT SPECIFICATION

CATALOG NUMBER: 1-028

PRODUCT DESCRIPTION:	Purified Human Creatine Kinase, MB Isoenzyme
SOURCE:	Human heart tissue
FORM:	Liquid
PURITY:	>90% by 10-20% gradient SDS PAGE 1 CK band by Corning agarose gel isoenzyme electrophoresis
BUFFER:	10 mM Tris, 5 mM mercaptoethanol, 50% (v/v) glycerol, pH 7.2 (at 25°C ± 3°C)
PRESERVATIVE:	None
STORAGE TEMPERATURE:	-20°C or colder
HAZARD/BIOHAZARD:	There is no known hazard associated with this product. Donor's serum was tested and found negative for HBsAg, HIV-1, HIV-2 and HCV antibodies by FDA approved methods.