



PRODUCT SPECIFICATION

CATALOG NUMBER: 1-029-ACT

PRODUCT:	PSA-ACT Complex
SOURCE:	PSA: Human seminal fluid (Donor serum tested and found negative for HbsAg and approved tests). ACT: Human plasma (Donor plasma tested and found negative for HbsAg and antibodies to HIV-1, HIV-2, and HCV by FDA approved tests).
FORM:	Liquid
PURITY:	>95% by SDS-PAGE
BUFFER:	10 mM sodium acetate pH 5.6, 0.15 M sodium chloride
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C short term (up to 3 weeks) -20 degrees C or colder long term