



PRODUCT SPECIFICATION

CATALOG NUMBER: 1-030

PRODUCT:	Purified Human Skeletal Troponin, T subunit
SOURCE:	Human pectoral muscle Donors' serum tested and found negative for HBsAg, and HIV antibodies by an FDA approved test.
FORM:	Liquid
PURITY:	>90% by SDS-PAGE
BUFFER:	10 mM phosphate, 0.5 M potassium chloride, 0.15 M sodium chloride, pH 7.2
PRESERVATIVE:	None
STORAGE:	2-8 degrees C short term (up to one month) -20 degrees C or colder long term
NOTE:	Troponin T is only sparingly soluble at neutral pH and physiological salt concentrations.