



## PRODUCT SPECIFICATION

CATALOG NUMBER: A05013

PRODUCT DESCRIPTION:	NGAL Monoclonal Antibody
IMMUNOGEN:	Human Neutrophil gelatinase-associated lipocalin
SOURCE:	Mouse Ascites
PURIFICATION:	DEAE chromatography
PURITY:	>95% by HPLC and SDS-PAGE
BUFFER:	0.015 M Potassium Phosphate buffer containing 0.85% NaCl, pH 7.2
PRESERVATIVE:	0.05% sodium azide
SUBCLASS:	IgG <sub>1</sub>
APPLICATIONS:	EIA (This antibody can be used as a detection antibody in a Human NGAL ELISA in combination with the monoclonal antibody A05012, A05014 and A05015).
STORAGE:	2-8°C short term; -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.