



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>A10461</b>   |
| product description: | AMH (Anti-Mullerian Hormone) Monoclonal Antibody, Purified  |
| lot number:          | [lot specific]  |
| immunogen:           | Human AMH   |
| source:              | Murine ascites  |
| form:                | Liquid  |
| purification method: | DEAE chromatography   |
| purity:              | >95% by HPLC and SDS-PAGE   |
| buffer:              | 0.015M potassium phosphate, 0.85% NaCl, pH 7.2  |
| preservative:        | 0.05% NaN <sub>3</sub>  |
| isotype:             | IgG <sub>2a</sub>   |
| concentration:       | [lot specific] mg/mL (OD280 nm, E <sup>0.1%</sup> = 1.4)  |
| specificity:         | This antibody reacts with human AMH.  |
| pairing:             | This antibody can be used in quantitative EIA and paired with A10460503P or A10462503P.   |
| storage:             | 2-8°C (one year); -20°C long term avoiding multiple freeze thaw cycles.   |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |