



PRODUCT SPECIFICATION

CATALOG NUMBER: A11021

PRODUCT:	Monoclonal Antibody to AFP
SOURCE:	Mouse ascites
IMMUNOGEN:	Human AFP
FORM:	Liquid in 15 mM Phosphate buffer Containing 0.15M NaCl
PURIFICATION:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
PRESERVATIVE:	0.09% Sodium Azide pH 7.2
SUBCLASS:	IgG2a
STORAGE CONDITIONS:	2-8°, Maintain -70°C long term. Avoid freeze /thaw cycles
APPLICATIONS:	Sandwich ELISA (This antibody can be used as a capture antibody in a AFP ELISA combination with the monoclonal antibody (A11022)
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.