



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A11070**

|                     |  |
|---------------------|--|
| PRODUCT:            | Monoclonal Antibody to Alpha Fetoprotein (AFP)   |
| IMMUNOGEN:          | AFP purified from human fetal cord serum   |
| FORM:               | Chromatographically purified in 0.05 M PBS, 0.1% NaN <sub>3</sub> ; pH 7.4   |
| AFFINITY CONSTANT:  | 6 X 10 <sup>9</sup>  |
| SUBCLASS:           | IgG  |
| STORAGE CONDITIONS: | 2-8°C  |
| HAZARD/BIOHAZARD:   | This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |