



PRODUCT SPECIFICATION

CATALOG NUMBER: A15090

| | |
|---------------------|---|
| PRODUCT: | Apolipoprotein A-1 Monoclonal Antibody |
| IMMUNOGEN: | Human Apo A-1 |
| PURITY: | >95% by SDS-Page |
| FORM: | Purified IgG from mouse ascites |
| BUFFER: | Phosphate buffer saline, pH 7.0 |
| PRESERVATIVE: | 0.02% NaN ₃ |
| SUBCLASS: | IgG ₁ |
| STORAGE CONDITIONS: | 2-8°C short term; -20°C long term |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product. |