



PRODUCT SPECIFICATION

CATALOG NUMBER: A15100

PRODUCT:	Apolipoprotein A-1 Monoclonal Antibody
IMMUNOGEN:	Human Apo A-1
PURITY:	>95% by SDS-Page
FORM:	Purified IgG from mouse ascites
BUFFER:	Phosphate buffer saline, pH 7.0
PRESERVATIVE:	0.02% NaN ₃
SUBCLASS:	IgG _{2a}
STORAGE CONDITIONS:	2-8°C short term; -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.