



PRODUCT SPECIFICATION

CATALOG NUMBER: A15150

PRODUCT:	Apolipoprotein B Monoclonal Antibody
IMMUNOGEN:	Human Plasma
FORM:	Immunoglobulin purified from mouse ascites in 0.01M tris, 0.05M NaCl, pH 7.2 with 0.05% NaN ₃ as preservative.
PURIFICATION:	Ion Exchange Chromatography
SUBCLASS:	IgG _{2b}
AFFINITY CONSTANT:	1 x 10 ⁸
STORAGE CONDITIONS:	4°C short term, -20°C long-term
CROSS REACTIVITY:	No cross reactivity with Apo A-I or Apo A-II