



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A15160**

PRODUCT:	Monoclonal Antibody to Apolipoprotein B
IMMUNOGEN:	Human Plasma
FORM:	Immunoglobulin purified from mouse ascites
PURIFICATION:	Protein A
BUFFER:	0.05 PBS, with 0.05% NaN <sub>3</sub> as preservative, pH 7.2
AFFINITY CONSTANT:	1.0 X 10 <sup>8</sup>
SUBCLASS:	IgG <sub>2b1</sub>
STORAGE CONDITIONS:	4°C short term, -20°C long term
CROSS REACTIVITY:	No cross reactivity with Apo A-1 or Apo A-II.