



PRODUCT SPECIFICATION

CATALOG NUMBER: A15170

PRODUCT:	Monoclonal Antibody to Apolipoprotein B
IMMUNOGEN:	Human Plasma
FORM:	Immunoglobulin from purified mouse ascites
PURIFICATION:	Protein A
BUFFER:	0.01M Tris HCl, 0.5M NaCl, pH 8, with 0.1% NaN ₃ as preservative
AFFINITY CONSTANT:	1 x 10 ⁸
SUBCLASS:	IgG _{2b}
STORAGE CONDITIONS:	4°C short term, -20°C long term
CROSS REACTIVITY:	No cross reactivity with Apo A-I or Apo A-II.