



PRODUCT SPECIFICATION

CATALOG NUMBER: A15190

PRODUCT:	Monoclonal Antibody to Apolipoprotein B
IMMUNOGEN:	Human Plasma
FORM:	Immunoglobulin purified from mouse ascites
BUFFER:	DEAE-chromatography purified IgG in 0.015 M potassium phosphate buffer, 0.15 M NaCl, pH 7.2, with 0.1% NaN ₃ as preservative.
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.