



PRODUCT SPECIFICATION

CATALOG NUMBER: A15370

PRODUCT:	Albumin Monoclonal Antibody
FORM:	Each vial contains antibodies in PBS, pH 7.4 with 0.1% sodium azide as preservative.
ISOTYPE:	IgG ₁
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity tested by electrophoresis.
STORAGE TEMPERATURE:	2-8°C
SPECIFICITY:	Human albumin. No cross-reactivity with other human proteins tested. Does not react with bovine serum albumin and egg white albumin.
APPLICATION:	ELISA
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.