



PRODUCT SPECIFICATION

CATALOG NUMBER: A26040

PRODUCT:	Carcinoembryonic Antigen (CEA) Monoclonal Antibody
IMMUNOGEN:	Human Tumor
FORM:	Protein A purified Ig from mouse ascites in 0.01M Na Phosphate, 0.15M NaCl, 0.05% NaN ₃ , pH 7.20
TITER:	2 x 10 ¹⁰ Assays/mg (RIA)
AFFINITY CONSTANT:	8 x 10 ⁹
SUBCLASS:	IgG ₁
CROSS REACTIVITY:	<1% NCA
PART OF A COOPERATIVE PAIR:	Yes, with clone A26050
REFERENCE ASSAY:	RIA
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.