



PRODUCT SPECIFICATION

catalog number:	A33130
product description:	Human Ventricular Slow Skeletal Myosin Light Chain 1 (MLC-1) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified ventricular/slow MLC-1
source:	Murine; bioreactor cell culture
form:	Frozen liquid
purification method:	Protein G chromatography
buffer:	PBS, pH 7.3
isotype:	IgG _{1, kappa}
concentration:	[lot specific] mg/mL
pairing;	This antibody can pair in immunoassays with A33150.
expiration date:	[lot specific]
storage:	-20°C avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.