



PRODUCT SPECIFICATION

CATALOG NUMBER: A42500

PRODUCT:	Barbiturate Monoclonal Antibody
PURIFICATION:	Protein A affinity chromatography
PURITY:	>98% by SDS PAGE
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.05% Sodium Azide
SUBCLASS:	IgG _{1,κ}
STORAGE CONDITIONS:	2-8°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.