



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A45540**

PRODUCT:	Alpha 2 Macroglobulin Monoclonal Antibody
IMMUNOGEN:	Purified human Alpha 2 Macroglobulin
FORM:	DEAE column chromatography purified immunoglobulin
BUFFER:	0.015 M KPO <sub>4</sub> buffer, 0.85% NaCl pH 7.40
PRESERVATIVE:	0.1% sodium azide
AFFINITY CONSTANT:	N/A
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITION:	2-8° C. DO NOT FREEZE.
PAIRING INFORMATION:	This antibody pairs with rabbit anti-PSA (#D63210259) in PSA-Alpha 2 Macroglobulin complex immunoassays