



PRODUCT SPECIFICATION

CATALOG NUMBER: A45570

PRODUCT:	Alpha 2 Macroglobulin Monoclonal Antibody
IMMUNOGEN:	Human Alpha 2 Macroglobulin
FORM:	DEAE column chromatography purified immunoglobulin
BUFFER:	0.015 M KPO ₄ buffer, 0.85% NaCl, pH 7.40
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁
AFFINITY CONSTANT:	N/A
STORAGE CONDITION:	2-8° C. DO NOT FREEZE.
PAIRING INFORMATION:	This antibody pairs with rabbit anti-PSA (#D63210) in PSA-Alpha 2 Macroglobulin complex immunoassays.