



PRODUCT SPECIFICATION

CATALOG NUMBER: A45590

PRODUCT:	Human Alpha-2 Macroglobulin Monoclonal Antibody
IMMUNOGEN:	Purified Human Alpha-2 Macroglobulin
FORM:	DEAE chromatography purified
BUFFER:	0.015M potassium phosphate buffer, 0.15M NaCl, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG _{2b}
STORAGE CONDITIONS:	2-8°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.