



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A46300**

PRODUCT:	Anti-CA 19-9 Monoclonal Antibody
SOURCE:	BALB/C Mouse ascites
IMMUNOGEN:	purified CA 19-9 from tumor cells
FORM:	Liquid
PURIFICATION:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
BUFFER:	15 mM Phosphate buffer containing 0.15M NaCl, pH 7.2
PRESERVATIVE:	0.09% Sodium Azide
SUBCLASS:	IgG1k
STORAGE CONDITIONS:	2-8°, Maintain -70°C long term. Avoid freeze-thaw cycles.
APPLICATIONS:	Sandwich ELISA (This antibody can be used as Bboth capture and label)