



PRODUCT SPECIFICATION

catalog number:	A50055
product description:	MDPV Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purity:	>90% by SDS-PAGE
buffer:	PBS, pH 7.4
preservative:	0.02% NaN ₃
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for methylenedioxypropylvalerone (MDPV).
pairing:	This antibody can be paired with MDPV-BSA V50050.
notes/applications:	This monoclonal is known to work well in lateral flow.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.