



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>A51820</b>   |
| product description: | Clonazepam Monoclonal Antibody, Purified  |
| lot number:          | [lot specific]  |
| immunogen:           | 7-amino-clonazepam-KLH  |
| source:              | Murine  |
| form:                | Liquid  |
| purification method: | Protein A chromatography  |
| buffer:              | PBS, pH 7.4   |
| preservative:        | 0.1% NaN <sub>3</sub>   |
| concentration:       | [lot specific] mg/mL  |
| pairing:             | This antibody can be paired with V51800, clonazepam-BSA.  |
| notes/applications:  | This monoclonal is known to work well for lateral flow.   |
| storage:             | 2-8°C short term (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.   |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |