



PRODUCT SPECIFICATION

CATALOG NUMBER: A51900

PRODUCT:	Benzodiazepine Monoclonal Antibody
IMMUNOGEN:	Benzodiazepine-BTG
FORM:	Protein A purified immunoglobulin from mouse ascites
BUFFER:	10 Mm sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% sodium azide
SUBCLASS:	IgG ₁ , kappa
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Cross reactivity was determined by ELISA with competition assays. Reacts with Benzodiazepine compounds and metabolites
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.