



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A52015**

PRODUCT:	Benzoylcegonine Monoclonal Antibody
IMMUNOGEN:	Benzoylcegonine-BTG
FORM:	Protein A purified immunoglobulin from mouse ascites
PURITY:	≥95%
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% Sodium Azide
SUBCLASS:	IgG1
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.