



## PRODUCT SPECIFICATIONS

**CATALOG NUMBER: A53132**

PRODUCT:	Amphetamine Monoclonal Antibody
IMMUNOGEN:	Amphetamine-BSA
FORM:	Protein A purified immunoglobulin from mouse ascites
BUFFER:	PBS
PRESERVATIVE:	0.02% Sodium Azide
ISOTYPE:	IgG <sub>1k</sub>
STORAGE CONDITIONS:	2-8°C
CROSS REACTIVITY:	Qualified for use in lateral flow applications. No cross reactivity with Methamphetamine, MDMA, Ephedrine or Pseudoephedrine
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.