



PRODUCT SPECIFICATION

CATALOG NUMBER: A62020

PRODUCT:	Adenovirus Monoclonal Antibody
IMMUNOGEN:	Adeno 3
FORM:	Protein A purified immunoglobulin from ascites
BUFFER:	0.01M Phosphate buffer, 0.15M Sodium Chloride, pH 7.1
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG _{1,k}
STORAGE CONDITIONS:	2-8°C short term -20°C long term
SPECIFICITY:	Specific for all 41 serotypes of Adenovirus
PAIRING INFO:	The A62020 is the detection antibody and pairs with the A62050 the capture antibody.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.