



PRODUCT SPECIFICATION

CATALOG NUMBER: A63100

| | | |
|---------------------|--|--------------------|
| PRODUCT: | <u>Candida albicans</u> Monoclonal Antibody | |
| IMMUNOGEN: | <u>C. albicans</u> Strain 36802 | |
| FORM: | Immunoglobulin purified from mouse ascites | |
| PURIFICATION: | Protein A | |
| BUFFER: | 0.02M phosphate buffer, 0.25M sodium chloride, pH 7.6 filtered with a 0.45um filter. | |
| PRESERVATIVE: | 0.1% sodium azide | |
| SUBCLASS: | IgG _{2a} | |
| STORAGE CONDITIONS: | 2-8°C | |
| SPECIFICITY: | <u>Strain</u> | <u>ELISA Titer</u> |
| | 36802 | 1:200,000 |
| | 26310 | 1:100,000 |
| | BCH054 | 1:25,000 |
| | SSH022 | 1:100,000 |
| | Negative to Trichimonas | |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. | |