



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A63130**

PRODUCT:	<u>Candida albicans</u> Monoclonal Antibody
IMMUNOGEN:	C. <u>albicans</u> Strain 36802
FORM:	Immunoglobulin purified from mouse ascites
PURIFICATION:	Protein A
BUFFER:	0.02M phosphate buffer, 0.25M sodium chloride
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Negative to <u>Trichomonas</u>
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.