



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A63520**

PRODUCT:	Rotavirus Monoclonal Antibody
IMMUNOGEN:	Purified Bovine Rotavirus cell culture adapted strain MR.
FORM:	Purified immunoglobulin from mouse ascites
PURIFICATION:	Protein A Sepharose >95% pure by electrophoresis
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>2a</sub>
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	P42 inner-capsid antigen