



PRODUCT SPECIFICATION

catalog number:	A63610
product description:	Malaria pLDH Monoclonal Antibody, Purified
lot number:	[lot specific]
form:	Liquid
purification method:	Protein G chromatography
purity:	>95% by SDS-PAGE
buffer:	PBS, pH7.2
isotype:	IgG _{2b}
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific to pLDH from <i>P. falciparum</i> , <i>P. vivax</i> , <i>P. malariae</i> and <i>P. ovale</i> .
pairing:	This antibody can be used as the conjugate antibody and paired with A63620, A63630, A63640 or A63650 as the capture antibody.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.