



## PRODUCT SPECIFICATION

catalog number:	<b>A63620</b>
product description:	Malaria <i>P. falciparum</i> LDH Monoclonal Antibody, Purified
lot number:	[lot specific]
form:	Liquid
purification method:	Protein G chromatography
purity:	>95% by SDS-PAGE
buffer:	PBS, pH7.2
isotype:	IgG <sub>2b</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific to pLDH from <i>P. falciparum</i> .
pairing:	This antibody can be used as the capture antibody and paired with A63610 as the conjugate antibody.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.