



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A66000**

PRODUCT:	ACTH C-Terminal Monoclonal Antibody
SOURCE:	Mouse
IMMUNOGEN:	ACTH-BSA
PURIFICATION METHOD:	Protein A affinity chromatography.
BUFFER:	PBS
PRESERVATIVE:	10mM NaN <sub>3</sub>
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8°C short term -20°C long term
SPECIFICITY:	Specific for CLIP (17-39 ACTH). This antibody has no cross reaction with Synacthen. Positive reactivity with rat C-terminal ACTH.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.