



PRODUCT SPECIFICATION

CATALOG NUMBER: A66020

PRODUCT:	ACTH N-Terminal Monoclonal Antibody
SOURCE:	Tissue culture supernatant
FORM:	1 mg of monoclonal antibody (A_{280nm}) dissolved in 1 ml phosphate buffered saline. Contains 10 mM sodium azide as preservative and 1mg/ml bovine serum albumin.
PURIFICATION METHOD:	Protein A Sepharose affinity chromatography.
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Synacthen (1-24 ACTH) and no cross reaction with CLIP (ACTH 17-39)