



PRODUCT SPECIFICATION

CATALOG NUMBER: A78040

| | |
|---------------------|--|
| PRODUCT: | Beta-2-Microglobulin Monoclonal Antibody |
| IMMUNOGEN: | Human Beta-2-Microglobulin |
| FORM: | Purified IgG from mouse ascites |
| PURIFICATION: | Protein A chromatography |
| PURITY: | 98% by SDS-PAGE |
| BUFFER: | 50 mM Sodium Phosphate, 500 mM Sodium Chloride, 10 % glycerol (w/v), pH 8.0 |
| PRESERVATIVE: | 0.05% sodium azide |
| SUBCLASS: | IgG1, Kappa |
| STORAGE CONDITIONS: | 2-8°C, Do not freeze |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |