



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A78050**

PRODUCT:	B-2 Microglobulin Monoclonal Antibody
SOURCE:	Mouse ascites
IMMUNOGEN:	<u>Human B-2 Microglobulin</u>
FORM:	Liquid
BUFFER:	15 mM Phosphate buffer containing 0.15M NaCl
PRESERVATIVE:	0.09% Sodium Azide pH 7.2
SUBCLASS:	IgG1k
STORAGE CONDITIONS:	2-8°, Maintain -70°C long term. Avoid freeze-thaw cycles.
APPLICATIONS:	Sandwich ELISA (This antibody can be used as a capture antibody in a Human B-2 Microglobulin ELISA combination with the monoclonal antibody A78051489P)
PURIFICATION:	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)