



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A83840**

PRODUCT:	Thyroid Peroxidase Monoclonal Antibody
IMMUNOGEN:	Human thyroid peroxidase.
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity is tested by electrophoresis.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
ISOTYPE:	IgG <sub>1</sub>
SPECIFICITY:	Human thyroid peroxidase
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.