



PRODUCT SPECIFICATION

CATALOG NUMBER: A94092

PRODUCT:	<u>Helicobacter pylori</u> Monoclonal Antibody
IMMUNOGEN:	Purified <u>Helicobacter pylori</u> antigen
FORM:	Purified IgG from mouse ascites
PURITY:	>95% by SDS-Page
BUFFER:	Phosphate buffer saline, pH 7.0
PRESERVATIVE:	0.02% sodium azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C short term; -20°C long term
PAIRING INFORMATION:	This antibody forms a suitable pair when used as conjugate antibody with clone A94090 as capture.
HAZARD/BIOHAZARD:	This product contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.