



PRODUCT SPECIFICATION

CATALOG NUMBER: D-402-C

PRODUCT DESCRIPTION:	Donkey anti-goat IgG Fc Polyclonal, Affinity Purified
FORM:	Liquid
PURITY:	Affinity pure
BUFFER:	10 mM sodium phosphate, 0.85% NaCl (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE TEMPERATURE:	2-8°C
SPECIFICITY:	This antibody reacts with heavy chains on goat IgG but not with light chains. No reactivity is observed to non-immunoglobulin goat serum proteins.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal. There is no known biohazard associated with this