



PRODUCT SPECIFICATIONS

CATALOG NUMBER: D-501-C

PRODUCT:	Donkey anti-sheep IgG (H & L), Affinity Pure
APPEARANCE:	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8°C short term, -20°C or colder long term
SPECIFICITY:	This antibody reacts with heavy chains on sheep IgG, and with light chains on all sheep immunoglobulins. No reactivity is observed to non-immunoglobulin sheep serum proteins.
ANIMALS:	Serum donors were certified disease free by a registered veterinarian in the U.S.