



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D25010**

PRODUCT DESCRIPTION:	Rabbit anti Cortisol										
HOST ANIMAL:	Rabbit										
IMMUNOGEN:	Cortisol-3-BSA										
PRESERVATIVES:	None										
AFFINITY CONSTANT:	$9.0 \times 10^9$ L/M										
TITER:	15,000 test/ml in $^{125}\text{I}$ RIA										
STORAGE TEMPERATURE:	2-8°C short term; -20°C long term										
CROSS REACTIVITY:	<table><tr><td>Cortisol</td><td>100%</td></tr><tr><td>11- deoxycortisol</td><td>5.7%</td></tr><tr><td>Prednisolone</td><td>36%</td></tr><tr><td>Corticosterone</td><td>3.3%</td></tr><tr><td>Cortisone</td><td>&lt;0.7%</td></tr></table>	Cortisol	100%	11- deoxycortisol	5.7%	Prednisolone	36%	Corticosterone	3.3%	Cortisone	<0.7%
Cortisol	100%										
11- deoxycortisol	5.7%										
Prednisolone	36%										
Corticosterone	3.3%										
Cortisone	<0.7%										
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.										