



PRODUCT SPECIFICATION

CATALOG NUMBER: D41030

PRODUCT DESCRIPTION:	Goat anti Adenovirus Antigen
IMMUNOGEN:	Gradient purified, disrupted Adenovirus type 5 virions
FORM:	Neat antisera
PURIFICATION:	Purified IgG fraction
PURITY:	>95%
BUFFER:	0.01M phosphate saline buffer, pH 7.2
PRESERVATIVE:	0.1% sodium azide
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
HAZARD/BIOHAZARD:	Donor material tested negative and found negative for anti-HIV 1/2, anti-HCV and HBsAg antigen. There is no known hazard associated with this product.