



PRODUCT SPECIFICATION

catalog number:	D50042
product description:	K2 Polyclonal Antibody, Purified
lot number:	[lot specific]
source:	Rabbit
form:	Liquid
purification method:	Purified IgG fraction
buffer:	PBS, pH 7.1
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL
pairing:	This antibody can be paired with K2-BSA V50044.
notes/applications:	This polyclonal is known to work well for lateral flow.
storage:	2-8°C (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.