



PRODUCT SPECIFICATION

catalog number:	D51710
product description:	Zolpidem Polyclonal Antibody, Purified
lot number:	[lot specific]
source:	Rabbit
form:	Liquid
purification method:	Protein A purified IgG fraction
buffer:	0.05M phosphate, 0.85% NaCl, pH 7.15
preservative:	none
concentration:	[lot specific] mg/mL
pairing:	This antibody can be paired with V51710, zolpidem-HRP or V51720, zolpidem-BSA.
notes/applications:	This polyclonal is known to work well for lateral flow and ELISA.
storage:	2-8°C short term (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.