



PRODUCT SPECIFICATION

CATALOG NUMBER: D60630

PRODUCT:	Rabbit anti-Campylobacter jejuni
IMMUNOGEN:	ATCC strain 29428
PURITY:	Purified IgG fraction, purity >95%
BUFFER:	0.01M PBS, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
TITER:	>1:1,000 by ELISA
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term
SPECIFICITY:	Whole Cells
OTHER:	This antiserum has not been absorbed and may react with related microorganisms.