



PRODUCT SPECIFICATION

CATALOG NUMBER: D81220

PRODUCT DESCRIPTION:	Goat anti-Human C4 Serum
IMMUNOGEN:	Purified Human C4
BUFFER:	0.05M Tris, 0.5M NaCl, pH 7.5
PRESERVATIVE:	0.09% Sodium Azide
TOTAL PROTEIN:	7.6 g/dL
BECKER TITER:	5.8 mg/ml
STORAGE TEMPERATURE:	-20°C long term, 2-8°C short term
PROCESS:	Defibrinated, delipidized and adsorbed by solid phase chromatography as required.
SPECIFICITY:	Monospecific by IEP with neat antibody versus 2X pooled human serum and neat, pooled human plasma.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.