



PRODUCT SPECIFICATION

CATALOG NUMBER: D81610

PRODUCT DESCRIPTION:	Rabbit anti CRP, IgG Fraction
HOST ANIMAL:	Rabbit
ANTIGEN:	Purified human CRP
PURIFICATION:	DEAE Chromatographically
BUFFER:	0.015M KPO ₄ buffer, 100mM NaCl, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
STORAGE TEMPERATURE:	2-8°C, Do not freeze
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product