



PRODUCT SPECIFICATION

CATALOG NUMBER: D81620

PRODUCT DESCRIPTION:	Anti-Human CRP Serum
HOST ANIMAL:	Goat
ANTIGEN:	Purified human CRP
BUFFER:	20mM Tris, 250mM NaCl, pH 7.4
PRESERVATIVE:	15 mM NaN ₃
TOTAL PROTEIN:	Reference (by Biuret)
TITER:	Reference
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term
SPECIFICITY:	Non reactive when tested against normal human serum and single arc when tested against elevated CRP human serum by IEP.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.