



PRODUCT SPECIFICATION

CATALOG NUMBER: D82320

PRODUCT DESCRIPTION:	Anti-Human IgE Serum Polyclonal
HOST ANIMAL:	Goat
ANTIGEN:	Human IgE Serum
BUFFER:	0.05M Tris, 0.5M NaCl, pH 7.5
PRESERVATIVES:	0.1% NaN ₃
TOTAL PROTEIN:	2.9 g/dL
TITER:	0.50 mg/mL (Becker)
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term
SPECIFICITY:	No reactivity on IEP using undiluted antiserum tested against twice concentrated normal human serum. Single arc when tested against IgE elevated serum.