



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D83720**

PRODUCT DESCRIPTION:	Anti-Human Fibrinogen Polyclonal, IgG fraction
HOST ANIMAL:	Goat
ANTIGEN:	Human Fibrinogen from plasma
FORM:	Purified IgG fraction from pooled antisera
BUFFER:	Phosphate Buffered Saline
PRESERVATIVE:	0.1% Sodium Azide
TOTAL PROTEIN:	54 mg/ml
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
SPECIFICITY:	No reaction with normal human serum by immunoelectrophoresis or ouchterlony. Single line of identity by immunoelectrophoresis or ouchterlony gel against normal human plasma.
NOTE:	Antiserum has tested negative for blue tongue virus and tested negative. There is no known biohazard associated with this product.